

See Bottom of Page

Glass Fibre Filters

Normally used as a depth prefilter and placed directly on top of a membrane filter (please choose the prefilter diameter specified for the holder).

Type 13400 'standard' contains an acrylic latex binder and has a high loading capacity.

Type 13430 as 13400 but is thicker and more effective for 'dirty liquids'.

Type 13440 is a finer binder-free type recommended for the prefiltration of relatively clean solutions such as tissue culture media.



Type 13400 standard filters

F313/13	13mm	pack of 200
F313/25	25mm	pack of 200

Type 13430 extra thick standard filters

F314/13	13mm	pack of 200
F314/25	25mm	pack of 200


Type 13440 binder-free filters

F315/13	13mm	pack of 200
F315/25	25mm	pack of 200

Many other diameters available - -please ask

Forceps for Membrane Filters

Membrane filters should only be handled with suitable forceps to avoid contamination that can result from hand contact. These stainless steel forceps have smooth, flat tips. They are autoclavable and can be flamed.



T169 Stainless steel membrane forceps