

Please Scroll Down

Leica AC20 Ultrathin Section Stainer



The Leica EM AC20 automatic contrasting system for ultrathin sections ensures minimum user contact with reagents and reduced reagent consumption staining up to 20 grids per run. Use of a peristaltic pump and non-contact valves allow the reagents to travel directly through tubing to the grid chamber, resulting in fast, high-quality double staining and easy maintenance.

Safer more Efficient Working

Minimum user contact with reagents provides safer working conditions in the laboratory. 60 runs per one set of Ultrostains and short washing cycles result in faster, more efficient laboratory operation with associated cost savings.

- S534** Leica AC20 section stainer
- S534/1** Ultrostain 1 (4 x 200ml 0.5% uranyl acetate)
- S534/2** Ultrostain 2 (2 x 200ml 3% lead citrate)