

Please Scroll Down



## Leica EM IGL Immuno Gold Labelling System

The **automated** Leica EM IGL immunogold labeling system eliminates the need for tedious manual labelling of grids and improves the accuracy and reproducibility of results essential for TEM image quality. Reagents are automatically applied in the correct sequence to ultra-thin sections, using one of 99 possible preset programs.

Cross contamination from forceps and loops during labelling is minimised allowing the best possible result. High throughput labelling of **24 grids per run** saves time and allows users to concentrate on other laboratory responsibilities.

**G273** Leica EM IGL Immuno gold labelling system