
# 4 Filtration

#### See Bottom of Page

## Glass Syringe with Needle



3 ring syringes, glass/metal with needle. The three rings allow a good grip so the hand pressure is easier and less tiring. A viton O-ring ensures a leak-proof seal between glass and cylinder and chromium plated base. The plunger is made from type 1.4401 stainless steel. The syringes are supplied with two connectors, a male Luer lock for plastic syringe filter holders and a male Luer Linden connector for stainless steel filter holders. These syringes can be autoclaved at at 121°C for 30 minutes with the plunger removed.

S149	5ml glass syringe, 3 ring with needle	each
S150	10ml glass syringe, 3 ring with needle	each
S151	20ml glass syringe, 3 ring with needle	each
S152	50ml glass syringe, 3 ring with needle	each

#### All Glass Syringe

All glass syringe with **metal Luer fitting**, fully interchangeable, centre nozzle. This economical all glass syringe will handle the more aggressive solutions which tend to attack the piston seal of plastic syringes.

S206	2ml all glass syringe	pack of 10
S207	5ml all glass syringe	pack of 10
S208	10ml all glass syringe	pack of 10
S209	20ml all glass syringe	pack of 10

#### Plastic Syringe with Needle

Disposable syringe with male Luer lock connector. Can be used in combination with all syringe filter holders.

S143	1ml plastic syringe with needle	pack of 10
S144	2ml plastic syringe with needle	pack of 10
S145	5ml plastic syringe with needle	pack of 10
S146	10ml plastic syringe with needle	pack of 10
S147	20ml plastic syringe with needle	pack of 10
S148	50ml plastic syringe with needle	pack of 10



### Plastic Syringe

Plastic syringes made form polypropylene making them resistant to almost all solutions or mixtures used in electron microscopy.

S142	50ml plastic syringe	pack of 10
S060	20ml plastic syringe	pack of 25
S061	10ml plastic syringe	pack of 25
S062	5ml plastic syringe	pack of 25
S063	2ml plastic syringe	pack of 25
S141	1ml plastic syringe	pack of 25

## Disposable Syringe Needle

Fits all syringe filter holders where male Luer slip is required. **Sterile**, ready for immediate use and therefore especially suitable for use with disposable syringes or syringe filter holders to pierce rubber membranes etc. or to guide the filtrate into a container. Supplied in packs of 10

**S064**Coarse; 19G x 1" (2.54cm) 1.1mm nominal diameter x 40mm long**S065**Fine; 23G x 1" (2.54cm) 0.6mm nominal diameter x 25mm long

#### Stainless Steel Needle

Fits all filter holders and is autoclavable. To pierce rubber membranes or silicone caps or to guide the filtrate into a receiver e.g. a Petri dish

S257 Stainless steel needle for Luer slip



**)** 

4