


## Flexi-tip/Foil Tweezers



For some applications the bulk of a pair of normal tweezers, or even vacuum tweezers, is too great for accurate, gentle handling of some delicate devices, samples and live specimens. These tweezers are made from 0.12mm (0.035") anti-magnetic steel and are ideal for handling plant and animal specimens with positive control. Acid resistant. Length 100 mm, very fine pointed end

**T544** Flexi-tip/Foil Tweezer, sharp, straight tips

A broad tip is sometimes preferable for some applications. This tweezer has a rounded blunt Tip, 100mm length x 0.12mm (0.035") Thickness. Extremely low force stainless steel, anti-magnetic & acid resistant.

**T545** Flexi-tip/Foil tweezer, rounded end, 100mm long