SEM, TEM, STM Supplies

Specialised Hitachi SEM Stubs

S570	Hitachi 15mm Ø x 10mm M4 angled SEM sample stub 45° Al	each
S570/1	Hitachi 15mm Ø x 10mm M4 angled SEM sample stub 45° Al	pack of 10
S571	Hitachi 25mm Ø x 16mm M4 angled SEM sample stub 45° Al	each
S571/1	Hitachi 25mm Ø x 16mm M4 angled SEM sample stub 45° Al	pack of 10
S572	Hitachi 15mm Ø x 11mm M4 angled SEM sample stub 45° & 90° AI	each
S572/1	Hitachi 15mm Ø x 11mm M4 angled SEM sample stub 45° & 90° Al	pack of 10
S573	Hitachi 25mm Ø x 16mm M4 angled SEM sample stub 45° & 90° AI	each
S573/1	Hitachi 25mm Ø x 16mm M4 angled SEM sample stub 45° & 90° AI	pack of 10
S574	Hitachi 25mm Ø x 16mm M4 angled SEM stub, double 90° Al	each
S574/1	Hitachi 25mm Ø x 16mm M4 angled SEM stub, double 90° Al	pack of 10
S575	Hitachi 32mm Ø x 12mm M4 angled SEM sample stub 45° & 90° AI	each
S575/1	Hitachi 32mm Ø x 12mm M4 angled SEM sample stub 45° & 90° AI	pack of 10
S576	Hitachi 15mm Ø x 10mm M4 dish SEM sample stub, Al	pack of 5
S576/1	Hitachi 15mm Ø x 10mm M4 dish SEM sample stub, Al	pack of 10
S576/2	Hitachi 15mm Ø x 10mm M4 dish SEM sample stub, Al	pack of 50

S573 S575

Pin Stubs with SEM Sample Clips

The S-Clip sample holders comprise one or more spring loaded sample clips either on an SEM sample stub or holder. The spring loaded S-clips are specifically designed to hold thin specimens up to 2mm. Optional elevating (blocks) pillars are available to accommodate thicker samples. Washers can also be used to increase the clamping height.

The S-clips are made from conductive spring grade copper alloy mounted with a small M2 screw which allows the clip to rotate. Samples are clamped securely without the need for adhesives or conductive pastes. Sample changing is quick, easy and clean; contamination issues due to adhesives or glues are eliminated. They are ideal to hold Si chips and other thin samples. The S-Clip sample holders are available to hold samples horizontal (flat), at 45° or at 90°. The 45° and 90° are useful for holding cross sections of thin samples and wafer pieces. S-Clip sample holders are a cost effective alternative for the 1"-3" wafer holders. They can hold the smaller 1" (25.4mm), 2" (51mm) and 3"(76mm) wafers if it is allowed to hold the wafer down with clips:

38m stub diameter is suitable for 1" (25.4mm) wafers: 63m stub diameter is suitable for 2" (51mm) wafers and 100mm stub diameter is suitable for 3" (76mm) wafers. The S-Clips are also separately available to enable customers to make their own custom

made sample clip holders. Dimensions are 12.8mm long, tip is 1.6mm wide, height is 1.5mm with a 2.2mm Ø hole for mounting with an M2 screw. Made from 0.25mm thin spring grade copper alloy. Included with the S-clips are the M2 (x3) mounting screws.

S587 19mm Ø standard pin stub with 1 x S-Clip for FEI/Philips, LEO etc

S588 19mm Ø short (6mm) pin stub with 1 x S-Clip for Zeiss

S577/1 25mm Ø x 9mm pin stub with 1 clip for FEI/Philips, ZEISS/LEO, Cambridge, each S577/2 25mm Ø x 9mm pin stub with 2 clips for FEI/Philips, ZEISS/LEO, Cambridge, each S577/3 25mm Ø x 9mm pin stub with 3 clips for FEI/Philips, ZEISS/LEO, Cambridge, each

Shorter version

S578/1 25mm Ø x 6mm pin stub with 1 clip for ZEISS/LEO, SEM, FIBs S578/2 25mm Ø x 6mm pin stub with 2 clips for ZEISS/LEO, SEM, FIBs S578/3 25mm Ø x 6mm pin stub with 3 clips for ZEISS/LEO, SEM, FIBs

S579/1 32mm Ø x 9mm pin stub with 1 clip for FEI/Philips, ZEISS/LEO, Cambridge, each \$579/2 32mm Ø x 9mm pin stub with 2 clips for FEI/Philips, ZEISS/LEO, Cambridge, each S579/3 32mm Ø x 9mm pin stub with 3 clips for FEI/Philips, ZEISS/LEO, Cambridge, each

S580/1 38mm Ø x 9mm pin stub with 1 clip for FEI/Philips, ZEISS/LEO, Cambridge, each S580/2 38mm Ø x 9mm pin stub with 2 clips for FEI/Philips, ZEISS/LEO, Cambridge, each S580/3 38mm Ø x 9mm pin stub with 3 clip for FEI/Philips, ZEISS/LEO, Cambridge, each



each

each

each







10

SEM, TEM, STM Supplies



Pin Stubs with Sample Clips cont...

S581/1 50mm Ø x std pin stub with 1 clip for FEI/Philips, ZEISS/LEO, Cambridge each S581/2 50mm Ø x std pin stub with 2 clips for FEI/Philips, ZEISS/LEO, Cambridge each S581/3 50mm Ø x std pin stub with 3 clips for FEI/Philips, ZEISS/LEO, Cambridge each S581/4 50mm Ø x std pin stub with 4 clips for FEI/Philips, ZEISS/LEO, Cambridge each S581/8 50mm Ø x std pin stub with 8 clips for FEI/Philips, ZEISS/LEO, Cambridge each

S582/2 63mm Ø x std pin stub with 2 clips for FEI/Philips, ZEISS/LEO, Cambridge each S582/3 63mm Ø x std pin stub with 3 clips for FEI/Philips, ZEISS/LEO, Cambridge each S582/4 63mm Ø x std pin stub with 4 clips for FEI/Philips, ZEISS/LEO, Cambridge each

S583/1 100mm Ø x std pin stub with 1 clip for FEI/Philips, ZEISS/LEO, Cambridge each S583/2 100mm Ø x std pin stub with 2 clips for FEI/Philips, ZEISS/LEO, Cambridge each S583/3 100mm Ø x std pin stub with 3 clips for FEI/Philips, ZEISS/LEO, Cambridge each S583/4 100mm Ø x std pin stub with 4 clips for FEI/Philips, ZEISS/LEO, Cambridge each S583/8 100mm Ø x std pin stub with 8 clips for FEI/Philips, ZEISS/LEO, Cambridge each

S-Clips 45° & 90° on 25.4mm Ø & 12.7mm Ø Pin Stubs

S589 S-Clip sample holder with I clip 45° on 25.4mm Ø pin stub each **S590** S-Clip sample holder swivel mount with 1 S-Clip on 25.4mm Ø pin stub each S591 S-Clip sample holder with 2x S-Clips 45°/90° on 25.4 m Ø pin stub S592 S-Clip sample holder with 2x S-Clips double 90° on 25.4mm Ø pin stub **S593** S-Clip sample holder with 1x S-Clip at 90° on 12.7mm Ø pin stub **S594** S-Clip sample holder with 4x S-Clips at 90° on 12.7mm Ø pin stub **S595** S-Clip sample holder with 4x S-Clips at 45° & 90° on pin stub each S596 S-Clip sample holder with 10x S-Clips at 90° on 10 x 50mm pin stub each

each each each

each

S-Clips on M4 Hitachi Stubs

S-Clip sample holder with 1x S-Clip on 15mm Ø x 10mm stub, M4 each S-Clip sample holder with 1x S-Clip on 25mm Ø x 10mm stub, M4 each S598/2 S-Clip sample holder with 2x S-Clip on 25mm Ø x 10mm stub, M4 each S598/3 S-Clip sample holder with 3x S-Clip on 25mm Ø x 10mm stub, M4 each

S599 S-Clip sample holder with 1x S-Clip on 32mm Ø x 10mm stub, M4 each S599/2 S-Clip sample holder with 2x S-Clip on 32mm Ø x 10mm stub, M4 each S599/3 S-Clip sample holder with 3x S-Clip on 32mm Ø x 10mm stub, M4 each

S630 S-Clip sample holder with 1x S-Clip on 50mm Ø x 6mm stub, M4 each S630/1 S-Clip sample holder with 2x S-Clip on 50mm Ø x 6mm stub, M4 each S630/3 S-Clip sample holder with 3x S-Clip on 50mm Ø x 6mm stub, M4 each S630/4 S-Clip sample holder with 4x S-Clip on 50mm Ø x 6mm stub, M4 each

S631 S-Clip sample holder with 1x S-Clip on 63mm Ø x 6mm stub, M4 each S631/2 S-Clip sample holder with 2x S-Clip on 63mm Ø x 6mm stub, M4 each S631/3 S-Clip sample holder with 3x S-Clip on 63mm Ø x 6mm stub, M4 each S631/4 S-Clip sample holder with 4x S-Clip on 63mm Ø x 6mm stub, M4 each S631/8 S-Clip sample holder with 8x S-Clip on 63mm Ø x 6mm stub, M4 each

S632 S-Clip sample holder with 1x S-Clip on 100mm Ø x 6mm stub, M4 each S632/2 S-Clip sample holder with 2x S-Clip on 100mm Ø x 6mm stub, M4 each S632/3 S-Clip sample holder with 3x S-Clip on 100mm Ø x 6mm stub, M4 each S632/4 S-Clip sample holder with 4x S-Clip on 100mm Ø x 6mm stub, M4 each S632/8 S-Clip sample holder with 8x S-Clip on 100mm Ø x 6mm stub, M4 each

SEM, TEM, STM Supplies

10

S-Clips on M4 Hitachi Stubs 45° & 90°

S633 S-Clip sample holder 45° with 1x S-Clip on 25mm Ø stub M4eachS634 S-Clip sample holder 45°/90° with 2x S-Clip on 25mm Ø stub M4eachS635 S-Clip sample holder 90° with 4x S-Clip on 12.7x12.7mm Ø stub M4eachS636 S-Clip sample holder with 4x S-Clips at 45° & 4xS-Clips at 90°, M4eachS637 S-Clip sample holder with 10x S-Clips at 90°, 10x50mm,stub M4each

Spare S-Clip Screws & Elevating Blocks

S638 5 Copper alloy S-Clips with 5 each brass M2 x 3mm screwssetS639 6 each S-Clip elevating pillars, 4mm height, M2 x3mm male/femaleset

S-Clips on Flat JEOL Specimen Stubs

S641 S-Clip sample holder with 1x S-Clip on $32mm \emptyset x 10mm$ JEOL SEM stubeachS641/2 S-Clip sample holder with 2x S-Clip on $32mm \emptyset x 10mm$ JEOL SEM stubeachS641/3 S-Clip sample holder with 3x S-Clip on $32mm \emptyset x 10mm$ JEOL SEM stubeach

S642 S-Clip sample holder with 1x S-Clip on 50mm Ø x 10mm JEOL SEM stubeachS642/2 S-Clip sample holder with 2x S-Clip on 50mm Ø x 10mm JEOL SEM stubeachS642/3 S-Clip sample holder with 3x S-Clip on 50mm Ø x 10mm JEOL SEM stubeachS642/4 S-Clip sample holder with 4x S-Clip on 50mm Ø x 10mm JEOL SEM stubeach

S-Clips on 45° & 90° JEOL Stubs

S643 S-Clip sample holder with 1x S-Clip on 45° x 25mm Ø JEOL stubeachS644 S-Clip sample holder with 2x S-Clips at 45° & 90° on 25mm Ø JEOL stubeachS645 S-Clip sample holder with 2x S-Clips double 90° on 25mm Ø JEOL stubeach





The EM-Tec SEM stub tweezers securely grip SEM stubs and mounts to avoid touching them by hand to eliminate contamination issues. They are made from *anti-magnetic stainless steel* and are clean room compatible. They are excellent quality at an affordable price **and are available for:**

All grooved pin stubs with 12.7, 19, 25.4 and 32mm head diameters (45° handle angle) T600 anti-magnetic S/S tweezers for grooved pin stubs 12.7mm Ø 120mm long each T601 anti-magnetic S/S tweezers for grooved pin stubs 19mm Ø 120mm long each T602 anti-magnetic S/S tweezers for grooved pin stubs 25.4mm Ø 120mm long each T603 anti-magnetic S/S tweezers for grooved pin stubs 32mm Ø 115mm long each JEOL cylinder stubs with 9.5, 12.2, 25.4 and 32mm head diameters (30° handle angle) T604 anti-mag S/S grippers for 9.5mm Ø JEOL cylinder stubs 114mm long each T605 anti-mag S/S grippers for 12.2mm Ø JEOL cylinder stubs 120mm long each **T606** anti-mag S/S grippers for 15mm Ø JEOL cylinder stubs 120mm long each 115mm long T607 anti-mag S/S grippers for 25mm Ø JEOL cylinder stubs each T608 anti-mag S/S grippers for 32mm Ø JEOL cylinder stubs 115mm long each Hitachi cylinder stubs with 15, 25 and 32mm head diameters (30° handle angle) T609 anti-mag S/S grippers for 15mm Ø Hitachi & ISI cylinder stubs 120mm long each T610 anti-mag S/S grippers for 25mm Ø Hitachi & ISI cylinder stubs 120mm long each T611 anti-mag S/S grippers for 32mm Ø Hitachi & ISI cylinder stubs 120mm long each

