


## Dissecting Awl

Stainless steel spike in a choice of two different lengths milled aluminium handles.  
Length of spike: 28mm.


<b>T358</b>	Dissecting awl - large	pair
<b>T359</b>	Dissecting awl - small	pair



## The Perfect Loop

Using the 'Perfect Loop', you can easily transfer thin sections from the microtome boats to grids with a minimum of creasing. The outside diameter of the loop (3mm) is the same as the grid. This assists the attraction of the loop to the grid when in water, which attach together through the surface tension of the water. Larger 7mm diameter loops are available for LM work (large sections).

<b>L138</b>	Perfect loop, EM set	each
<b>L141</b>	Replacement EM loop only	each
<b>L142</b>	Perfect loop, EM handle	each




## Dissecting Sets

### Basic Set

A 10 piece set comprising 2 solid scalpels, 2 needles, 1 seeker, forceps-blunt and sharp, scissors-fine and blunt and a fine brush all in a PVC roll. (Needles and seeker in wooden handles).


<b>D078</b>	Dissecting set - basic 10 piece	each
<b>D078/S</b>	Dissecting set - 10 piece, loose blades	each



### Advanced Set

A 14 piece set comprising 3 solid scalpels, 2 needles, 1 seeker, 1 section lifter, forceps-blunt, sharp and sharp curved, scissors - fine and blunt, 1 section razor and a fine brush all in a PVC roll. (Needles and seeker in stainless steel handles).

<b>D190</b>	Dissecting set - advanced 14 piece	each
<b>D190/S</b>	Dissecting set - 14 piece, loose blades	each



### Options:

The **S** sets have the solid scalpels replaced with Scalpel handles no.3 & 4 and 20 blades. If blade type is not specified at ordering, TAAB will supply 10 x No.10 and 10 x No.22a

## Dissecting Rolls

A choice of two dissection rolls both with 13 pockets and tape tie, one is manufactured in khaki canvas and the other in PVC with transparent pockets.

<b>T360</b>	Dissecting roll - canvas	each
<b>T361</b>	Dissecting roll - PVC	each

