5

## **Specimen Preparation**

## **Cryo EM Pucks & Grid Box Holders**

Keep your samples safely organized during storage and transportation. This puck/grid storage design is an advanced system for easily keeping track of your Cryo-EM specimen grids under cryogenic conditions.

The Cryo-EM Pucks feature 12 numbered wells for circular grid boxes, and a seethrough lid with notches for easy tweezer access for retrieving the grid holders. Each puck can be customized with a unique serial number, so that every sample is identified ready for insertion into the microscope. When using cryo grid boxes with flat lids, 24 can be stored per puck.

**Our Cryo-EM Pucks now feature a Lid Lock mechanism.** This allows the pucks to be easily and safely loaded into storage and shipping carriers (canes). Shelved shipping cane can hold up to 7 pucks and shelved storage cane up to 10 pucks.

C439 Cryo EM Pucks, standard 12 hole each
Cryo EM Pucks, standard 12 hole set of 7



Carriers can be used for the transport of Universal V-1 Pucks, ALS Pucks & Cryo-EM Pucks in dryshippers, such as the CX100. This set includes a shelved shipping cane, locking rod and hooked handle (shown right). Seven pucks (Uni-pucks, ALS pucks or Cryo EM) may be shipped or stored in this cane.

C441 Shelved storage cane holds 10 Cryo Pucks (HC34 compatible) each
 C442 Shelved shipping cane (ALS style) with locking rod and hooked handle
 C443 Protective case for storage and shipping canes



C443 Protective case for Storage and Shipping canes


## **Cryo-EM Starter Kit - Six Full Storage Canes**

Organised storage for Cryo-EM specimen grids under cryogenic conditions. The Cryo-EM Pucks feature 12 numbered wells, a see-through lid and easy-access tweezer notches for easy grid box retrieval. Each puck has a unique serial number to keep track of every sample. This kit contains six sets of pucks and canes for Cryo-EM. The pucks and tools can also be purchased individually.

The Cryo-EM Full Kit is a complete set of pucks and tools for Cryo-EM including:

- 60 x Cryo-EM Pucks
- 6 x Shelved Storage Cane (HC34 Compatible)
- 1 x Angled Cryo-Tongs for pucks

C444 Cryo EM starter kit—6 full storage canes and tools









## C444 Starter Kit

Six complete sets of pucks for Cryo-EM storage and transport under cryogenic conditions





Phone: +44 (0) 118 981 7775 Fax: +44 (0) 118 981 7881 E-mail: sales@taab.co.uk